

KRTAP11-1-keratin associated protein 11-1 Gene

Reference: PTXBC130555

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HACL-1; HACL1; KAP11.1; keratin-associated protein 11-1; high sulfur keratin-associated protein 11.1; keratin associated protein 11-1		

For research use only.

GENE DESCRIPTION	keratin associated protein 11-1
GENE ID	337880
ORIGIN SPECIES	Human
ACCESSIONS	BC130555
VECTOR	pDONR223
NCBI SYMBOL	KRTAP11-1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcttcaactgtccacaagaattgtcttccaggccattggaggacgtgcaattgtccagtggcccaagttaccacgacttcaccactgatgtgactgcctggcgccgcatctgtttgccagttcctccagactggctcttggtcctctggaccactgtcaagagacgtgtgtgag
cccactgttgcagcccaactgtttaccggcgaacttcatgtgtctccaaccttgcagggtgactgtctcgacaaactacgtgtatttccaacctgtctcaactacacagccggcgtcaccttgtctctagtggaatgtcagccctgggagcattctccagtgtctgccaaccagt
ggcggcatctctactgtctgccaaccagtgggaggagtctctactgtctgcccagcagcctgtgggtctccaggacgtatcagcagtcctgcgtgtccagctgccgaagaacctgctaa