

C17orf87-chromosome 17 open reading frame 87 Gene

Reference: PTXBC130556

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	transmembrane protein C17orf87; C17orf87; UNQ5783; SLP adapter and CSK-interacting membrane protein; DTFT5783; SLP65/SLP76, Csk-interacting membrane protein; SLP adaptor and CSK interacting membrane protein		

For research use only.

GENE DESCRIPTION	chromosome 17 open reading frame 87
GENE ID	388325
ORIGIN SPECIES	Human
ACCESSIONS	BC130556
VECTOR	pDONR223
NCBI SYMBOL	C17orf87
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggatacttcacagttcaggattccactgcaatgagctgggtgaggaataattctgcatcatcttagctgtggccatcattgtgtctctgtggcctgggcctcatcctgtactgtgtctgtaagtggcagcttagacgaggcaagaatgggaattgccaagccctgaaacacaagcaa
gtatgatgaagaaagatgtatgagaatgttcttaatgagtcgccagttcaattaccgcctctgccaccgaggaattggcctctctagaagactctccccacaggaagccccaagtcagccgcccgtacatactcactgttaataaagttaaaaaataagaagactgtttccatcccaagcta
cattgagcctgaagatgactatgacgatgttgaatccctgcaataactgaaaaagcatcatttga