

## KIAA1524-KIAA1524 Gene

Reference: PTXBC130564

|                |                                                                                                                        |                     |                |
|----------------|------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                    | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                   | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | CIP2A; p90; protein CIP2A; cancerous inhibitor of PP2A; cancerous inhibitor of protein phosphatase 2A; p90 autoantigen |                     |                |

**For research use only.**

|                         |                                                                                                  |
|-------------------------|--------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | KIAA1524                                                                                         |
| <b>GENE ID</b>          | 57650                                                                                            |
| <b>ORIGIN SPECIES</b>   | Human                                                                                            |
| <b>ACCESSIONS</b>       | BC130564                                                                                         |
| <b>VECTOR</b>           | pDONR223                                                                                         |
| <b>NCBI SYMBOL</b>      | KIAA1524                                                                                         |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                   |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon)

atggactcactgctgctgaagtcctgctcctgactgctcagctacaaagccgtgaagtcagagggcgaacccactcagcttttgcggcacttggaggttaattctggacagaaactcacacgaactattacataaatcagatattaaagtgatgctgagttgcttagagcta  
cttgaagaccccaacataaagtcctcactgacttaagttattatcgggttgcgtcactagcagtagacattgaaccagagattgtcttcagaatacatataatctgaatagtgctggcgggagtggtttgctggagcagccacactgattcgggttttgcagtgactcaactctacaga  
agttacataataatgcaaaatttctattctggtgccaatagatgaatattacgttctcgtatagatcacattcaactcttgaagatgagttaaaaatgccttgcctagattattggcacaactcttctgttcaaacgcacataaagacattgagtaattgaaatcttttatcg  
aactcttatcactctgttggcccatagtagtttaactgtggttgcctttcaattatccagtttgacattaaatgaagaggtgggggaaaaagctattccatgctcgaacattcatcagacttttcaactaatatttaatttctcaaaacggtgatggcacttaactagaagtattcagttg  
acctactgagtgatctccttaagaatcctaaaaattgctgattatctcaccagatagagcacttttctcatgcttcccaagattataggcttcttcaatggaaagatcctgattccttcaaaaggttttagaattacttctgcttctgtcagtgactcagctgcgccatagctcactcagatgattg  
ttgaacagctccactgagcagccactctgggaagccataactaattgttagaactcactgtggtcctactgcgctgttgaagccaactttggacggatcagaaaaactgttctgttttagcattggagttgttcaaggaatattggagatgcatagatgctgtaactgttctcggctgac  
gttttggacccttctgctcctacaactcctgatcaactcagttcacagaacaaaatctagatgaggttcttaacagaaaaaatgtgaaggttgcgaagccattggaagttttgtaactctctgtggagatgatacactaaaaatgcatttgcacaaactgtcaagtgtacca  
ctcttagacaacaattacatatggcaagattgacgtgggatttgaacaaaggttcagacttctgaattatgcaaaactgctgctgatgtaatttgaacacttctgatttgattaaacaaacttaaacattgttcttgggtatggaagtaagcttcaaaaactctcaggaccacgtttgatta  
ctccttggctttgctttaaactcagataatagagaacaagtagactctggactgagaataattttagaggtgctcctcactgccagattttctgctttagtacttggagaagtagtagcagcaacaatgcctatagacaacaggaaacagacataaccagaaaaatgccctggcaatcatc  
aaatcacagtttccaacatcaataaaggttttaactcctcattgaaagatggttctcctgattgaattggaagaattatagagaacacttcagcttggaaatggttggtaaggtatcagatttggatgtgagaatattcgacataatggatgtatatgaaatgaactatccacattagcttccaaag  
aaagcaggttacaagatcttttgaacaaaaagctcagccctgcacaggtgatagactgattgtcagctcgtctgctcaaaactcaagctgaacagaggcagcagacactgctagtattgttagagaagttgagagaaaaatgaagagcttagtgtgttctgaaggcgcagca  
agttgaatcagaaaagcgcagagtgatattgagcatctttcaacataataggaattgagctgtggttgaagaacatgaaatactgacaaaactctacatggaacttctcagagaatgaaagtactgaaaagaagaataaagattacagatcacatgtgattcttgaataaacaat  
tgagacagtgaaaaagttgaatgagtcactcaaggacaacaaatgaaaaagtttgcctaattatagagaagaagacagagaaaaagagtaacagatcagctagtagacagagaacataagctagcgaatttgcatacaaaaaacaaagtacaagaagaagaattaaacacttaca  
aaaggaagggaagataagggaagaacattgataccttgaagaagaattaaagcagaacagataaagaagaaggttgagcattaaaggcttctcctcagaggttcaaaaggcacaattagaaggtctgttggaaagagaagagtccttggtaactcagcaagaggaattgaa  
caaacactcccatagatagcaatgatccacagtttaagtggttggaaaaataaattccagaacgtgtgaattcagtatatag

More info: Please contact us for any question: [contact@proteogenix.fr](mailto:contact@proteogenix.fr)