

C20orf118-chromosome 20 open reading frame 118 Gene

Reference: PTXBC130646

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C20orf118; TLD domain-containing protein 2; TBC/LysM-associated domain-containing protein 2; TBC/LysM-associated domain containing 2		

For research use only.

GENE DESCRIPTION	chromosome 20 open reading frame 118
GENE ID	140711
ORIGIN SPECIES	Human
ACCESSIONS	BC130646
VECTOR	pDONR223
NCBI SYMBOL	C20orf118
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagaggcctccgtggcgttacactcggctgccagcagggtggaggacacccctgtctggggaggagggtaacgaagagggaaggaggaggcagctccagaccagctgctgctcctgaggatccacgggtgccagctgacagaagccagccaggtttgagtgccct
cagagattcggcagctcagctttcacttcccacaaagagtcacggccatccctggagctgtgtttctgcacgtcaaggaggggtttcagcctgcagagcctgtaccggcgatggagggtgcagcgggccaagtctgtgtgtctcaggaccaggacggcagatattggagcc
ttctcctctcggctatccgactcagcaaggcttctatgtgactggcagacattctctctctctcccccacagctgaaggctcttaagtggactgggaagcaactcttcttgtgaaggaggactggattcactgatgatgggcagtggcagtgccgggttgggctgtgtgtgagtgag
actgttccgggggaagctcccctgccgacctcaacaacgaggtgctggccggcaggagcagttctgcatccaggagctggaggcttggtctcagctga