



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

KRTAP10-9-keratin associated protein 10-9 Gene

Reference: PTXBC131613

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	KAP10.9; KAP18.9; KRTAP18-9; keratin-associated protein 10-9; high sulfur keratin-associated protein 10.9; keratin associated protein 18-9; keratin-associated protein 10.9; keratin-associated protein 18.9; keratin associated protein 10-9		

For research use only.

GENE DESCRIPTION	keratin associated protein 10-9
GENE ID	386676
ORIGIN SPECIES	Human
ACCESSIONS	BC131613
VECTOR	pDONR223
NCBI SYMBOL	KRTAP10-9
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr