

KRTAP13-4-keratin associated protein 13-4 Gene

Reference: PTXBC131691

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	KAP13.4; keratin-associated protein 13-4; keratin associated protein 13-4		
For research use only.			

GENE DESCRIPTION	keratin associated protein 13-4
GENE ID	284827
ORIGIN SPECIES	Human
ACCESSIONS	BC131691
VECTOR	pDONR223
NCBI SYMBOL	KRTAP13-4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctacaactgctgcttagaaacttctctcccgctccttgggggtacctgtactaccaggtctctaccagcagcctgtgtctacagcactgcccctgtctctccagcactgccagctgcgttctctctctacagggactgtcagaagacctgctgggagccaccagctgcc
agaaatcctgtaccgccccagacctccatcctctgtctccctgtcagacgactgtctgtggtctctaggtcttccgtccagcagctgtcgtccagggctatggatctaggtgtctactcgtctgggaatggatccagtggtctcagattcctgaaatggaggtgtgttttcttc
cctgagttacggatccagattctgtaccacaaactactggtctctggagcctggcagctcttctgttacagaccaatctgtggatctcgtctatcaattcacctgctaa