

FLJ38723-hypothetical FLJ38723 Gene

Reference: PTXBC131692

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized FLJ38723		
For research use only.			

GENE DESCRIPTION	hypothetical FLJ38723
GENE ID	255180
ORIGIN SPECIES	Human
ACCESSIONS	BC131692
VECTOR	pDONR223
NCBI SYMBOL	FLJ38723
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtgacgtgtcacacagcaataaggagctgaagcaggagaactcagcttggcagagcagcttcaagttgtactgattgacaaggcaggatgcagtgtgacttgaggagctgaaaaagaaactggaactcactgagctcacgcttcaacagttatctagctggtgtgaa
gcccctgacgtaaccagcagttacagcagcccacagacgaacgggcacagctggaagcacacctggggcaggtgatggagtgctgaagtatctgcaaatggagagagagcagtatgccgagtatctgcatggagagagtgccatgtgtggcagagatgcgggagatgtcag
agcaggtgcacacgttgaggaggagagagtgacagcatgagtcgagtagcaggagctggagacgatctggctgaactaaggaaccagctgaacctgccccggagcctccagcagggccctccgaagtggagcagaagctacaagcagaggctgagcacctgtggaagga
gctggagaatctggcaggacagcttcaagcccaagtggagaaaatgaggcttga