

LCE3E-late cornified envelope 3E Gene

Reference: PTXBC131716

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	LEP17; late cornified envelope protein 3E; late envelope protein 17; late cornified envelope 3E		

For research use only.

GENE DESCRIPTION	late cornified envelope 3E
GENE ID	353145
ORIGIN SPECIES	Human
ACCESSIONS	BC131716
VECTOR	pDONR223
NCBI SYMBOL	LCE3E
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcctgccagcagaaccagaagcagtgccaacccccaccagtgccctcaccgaagtgtccccaaagaaccagctacagtgctgcctccagctcctctggctgtgccccaaagctctggggctgtggccctagctccgagggcggtcttctgaaccaccacagcgcca
ccaccgatgccggcgccagaggtccaactcctgtgacaggggcagtggtcagcaaggcggggctctgctgctgccacgggtctggggctgtgctga