

### SYT14L-synaptotagmin XIV-like Gene

Reference: PTXBC132697

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SYT14L; CHR14SYT; Strep14; syt14r; yt14r; synaptotagmin-16; CTD-2277K2.1; chr14 synaptotagmin; synaptotagmin 14-like protein; synaptotagmin XIV-related protein; synaptotagmin XVI; synaptotagmin 16		

**For research use only.**

---

<b>GENE DESCRIPTION</b>	synaptotagmin XIV-like
<b>GENE ID</b>	401135
<b>ORIGIN SPECIES</b>	Human
<b>ACCESSIONS</b>	BC132697
<b>VECTOR</b>	pDONR223
<b>NCBI SYMBOL</b>	SYT14L
<b>LENGTH</b>	ProteoGenix bp
<b>SEQUENCE PRIMERS</b>	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

**ORF (Only for reference, product provided without STOP codon )**

atgtctttagtgaaagtacatccgatgtcagtcctcgaagacgggtcagttccagaattctcattagcctgtttataatgccaaaatgggagactatcagcagaagtataaaaggcagccacttaaaaaattgggcagcaaacagaccaccaatacatatgttaagtaactctacg  
gaaatccacggatcaagagatgcctcaaatgcaagatccatccgcagagggcggcgaatccagtatacaaggaaacttttttcaaaagtaccctatttcagctttctcatgtgacactcatgctgtctgtgtatacaaaaagcagcatgagaagaagatgataggctggattatttagg  
tctcaacagctctggagaagaactcaatcactggactgaaatgaaagagtcaaaaggacggcaagtacgtagatggcagcgctgtctagatcacgatga