

HMX2-H6 family homeobox 2 Gene

Reference: PTXBC132758

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	homeobox protein HMX2; H6L; Nkx5-2; homeo box (H6 family) 2; homeobox protein H6 family member 2; H6 family homeobox 2		

For research use only.

GENE DESCRIPTION	H6 family homeobox 2
GENE ID	3167
ORIGIN SPECIES	Human
ACCESSIONS	BC132758
VECTOR	pDONR223
NCBI SYMBOL	HMX2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggcagcaagaagatgcgggcaaggggtgctggcgccgtggtccagctcaccatccagtcacatcctggcgggggccctcggaggcaccgcggggagcccgctggctggccagccaggaagcgagcctgtccgtgtcctcggaggaggaggagccggacg
acggctggaagcgccccctgtctgccccagaccagcacggccctaaggagcaggccccaagcaccatcccccatccctttcttgcctgggtaccaccaagggcagcgaggctcggcccgccggcggttgagcgcacgcctttctctccttcgactcggactttaa
gaagagaaaagaggctcctgccccgcggtcgcctcgcgggtgctcagcgccgagggcgccggggacggcggtgctgagcgagcgccggcgccgccaagaagaagacgcgcaccgtcttttcgcgcagccaggtgtaccagctcagtcacacttcgacatgaagcgtacactgag
cagctcggagcgccgtgctcgcctcagcctgcagctcacggagaccaggttaagacttggttccagaaccgcccgaacaagtgggaagcggcagctctcggctgagctggaggcgccaacatggcgacgcgtcggcgagactctggtgagcatgccgtgtgtccgg
gacagttcgtctgcgtgcccgtgcccgtcgtcgtcgtcttccccgcgcgtctactaccgggaagcaacctctcgcccttacctctctacaacctatacaacaagctcgactactga