

RAB27A-RAB27A, member RAS oncogene family Gene

Reference: PTXBC132800

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RAB27A, member RAS oncogene family; GS2; HsT18676; RAB27; RAM; ras-related protein Rab-27A; GTP-binding protein Ram; mutant Ras-related protein Rab-27A; rab-27		

For research use only.

GENE DESCRIPTION	RAB27A, member RAS oncogene family
GENE ID	5873
ORIGIN SPECIES	Human
ACCESSIONS	BC132800
VECTOR	pDONR223
NCBI SYMBOL	RAB27A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctgatggagattatgattacctcatcaagtttttagcttgggagactctggtgtaggggaagaccagtgactttaccaatatacatggttaatttaactccaaatttatcacaacagtgggcattgatttcagggaagaaagagtggtgtacagagccagtggccggatggagccactg
gcagagccagagaaatccactgcagttatgggacacagcagggcagggagaggtttcgtagcttaacgacagcgttcttcagagatgctatgggttttctctactttttgatctgacaaatgagcaaatgttcctcaatgtcagaaactggataagccagctacagatgcataatgtgaa
aaccagatatagtgctgtgtggaacaagagtgatctggaggaccagagatagtgaaagaggaggaagccatagcactgcagagaaataggaatccctactttgaaactagtgctgccaatgggacaaacataagccaagcaattgagatgctcttgacctgataatgaagcga
atggacgggtgtgtggacaagtcctggattcctgaaggagtggtgcgatcaaatggtcatgcctctacggatcagtttaagtgaagaaaggagaaaggggcatgtgctgttga