

C3orf42-chromosome 3 open reading frame 42 Gene

Reference: PTXBC132802

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C3orf42; GHRL-AS2; GHRLOS2; NAG73; long intergenic non-protein coding RNA 852		

For research use only.

GENE DESCRIPTION	chromosome 3 open reading frame 42
GENE ID	84657
ORIGIN SPECIES	Human
ACCESSIONS	BC132802
VECTOR	pDONR223
NCBI SYMBOL	C3orf42
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaccgagtgccctccccgggaagaccagtgtcgttaactggaattaagttgtggaacataagagtcaggctcgtgtctgctgtaactggagtgagggaatgtttgggatccacctgggtatcagtaagggaacatggcgacggccagcctgccccatgtgaggcatttgcttgcagggccttcagtagtgagaaggaggaggagaatgcagctcttcagtgctgtctagcagtgcccaagagccgtgattgggcacctcacttaacgtctaactcagattcactcttgggcattcctgtcttctctccaatctga