

C6orf221-chromosome 6 open reading frame 221 Gene

Reference: PTXBC132844

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C6orf221; ECAT1; HYDM2; KHDC3-like protein; ES cell-associated transcript 1 protein; KH domain containing 3 like, subcortical maternal complex member		

For research use only.

GENE DESCRIPTION	chromosome 6 open reading frame 221
GENE ID	154288
ORIGIN SPECIES	Human
ACCESSIONS	BC132844
VECTOR	pDONR223
NCBI SYMBOL	C6orf221
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggacgtccagcggttccgacgtctgcaactgatgcagccaaagcaatgccagtggaggtgctcggtcacccctaaagcgttctctggttcactctgagttcctgaagaatccgaagtagttgccttgaggtttggctgggaaaagatctcgccggggcgag
aacgatcccgacgtccaggtatgtccaaatctgattcagctgaatcgattggaccctaacggcgaggtgagatcttggatttgggagggcctttaccaggaggacacaaatgaatgatcatgaacctggctgactatcaccgccagctccaggcgaaggctcaggaaggc
cctgcccaggatgtgccactcagaaggccgagaccagcgggttcaatagaagtcggggagggccggacgcagcgttcggtggaggtccgggagggccggaccagcgttcgggtggaagtcaggaggtcgggacacaggggttccgggtggaggtgcaggagccggga
cccagcagtcctccaggtgccaacaagtcggggaccagcgtatccccgaagctgcagcgaaggcagtgacccagcgggttcgcgaggtatgccgggacccagttactagattatga