

C15orf28-chromosome 15 open reading frame 28 Gene

Reference: PTXBC132845

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C15orf28; HsT18971; NCRNA00321; ANP32A intronic transcript 1 (non-protein coding)		

For research use only.

GENE DESCRIPTION	chromosome 15 open reading frame 28
GENE ID	80035
ORIGIN SPECIES	Human
ACCESSIONS	BC132845
VECTOR	pDONR223
NCBI SYMBOL	C15orf28
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagcctatctgcagaggttggcctgaagccataatgagcagttctgttcttttagcccttgtctagtttatggctctgttcagacatgcccttctgtcttggtaaaagtactgggtgtgtctttgcatltaggtgtgtgtctgtgagctctgggaaagccttctagaggagagccctccagggttgc
tctttgggtgggagtagtgggagcggaaatggaacgtgactgccctctaaagccagaagggtacttcattcaggcgaggcaggttggcatgtgtagtcactcaggagcttctcagttcagcgccgggtgttcacatcaggagtgagcagggttgatggagtgagaatgtgcc
tctga