



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

FLJ40292-hypothetical LOC643210 Gene

Reference: PTXBC132855

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

For research use only.

GENE DESCRIPTION	hypothetical LOC643210
GENE ID	643210
ORIGIN SPECIES	Human
ACCESSIONS	BC132855
VECTOR	pDONR223
NCBI SYMBOL	FLJ40292
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcgggcattattccaaagccttttgcgtaggggcgcctccagggtccacggcagcctgtgcagcatcggccgtcgcccaaggtgatacaggtgagggccgtgagagggggccaggtggcgccctccagactcactgcctctgtgtaggaaggcgaggaatgggagcctgttccaaaggtcagccgccacgcctgtgtcggggagcattggcgcaaacacggggagggagggcggcagagccctggcaaacagacactcagtttcgcacatggccctcaagccgtgtgcctgttggcgtgtggcgtgggggtttctcattcagcaggtctgcctgcctgctggattttgctttctatagatgttccccgtttctctcatcgtctctgcacccacacttttgcgaacacactgttgcctctctctctctctctctg