



Reference: PTXBC132898

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C13orf13; PRO28631; TMEM46; WGAR9166; bA398O19.2; hShisa; protein shisa-2 homolog; shisa homolog 2; transmembrane protein 46; shisa family member 2		

For research use only.

GENE DESCRIPTION	shisa homolog 2 (<i>Xenopus laevis</i>)
GENE ID	387914
ORIGIN SPECIES	Human
ACCESSIONS	BC132898
VECTOR	pDONR223
NCBI SYMBOL	SHISA2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

atgttggcgcctcgcgcctcgtcgtctctcctcctcgtcgaacccgcctcgtcgtctcgcagctcgtcgtcgtcgtcgtcgcgcgcggggcgagcgccgagccgcgcgcgagctgccaccgctgctcgtcggacgcgcagcgcgcgtctgcgcgcctcgcctccagtgtcccgagcgtctgcac
ggcggcgacgccaccatctcgtcgtcggcagctcgcggttgcgtactcgtcgtccagcggcggaggcgcgcgtcgtgaccagggcggtcgtcgacaatgaccgccagcaggcgctgtgcgagcctgcgccggcgccgacaaagcggccccgagcgtcgcgagtgcccatctacgtgcc
gttctcattgttgcgtccgtgtttcgtcgtcttaccatctgtggccctcgtggccgcctcgttgcgtcagatgtctcgcagcaacgagcatcccacagagcgcgcgcagccggcgagtcaccgtctgatggccacatcccattgatccccgagtcgcacgacactccccggcgtctcctc
acgcagcttcgcgcagacgctgcgcaggttcagctccagccgcaactcagggcgccccggcccccacaagctgcagacacaaactcgttgcgtcggcgaagggagacatgaacaacgtgtatgtctcaacatctgcgcagcaalttctcgtcgtgaactgcgcagccaccagattgtgcac
atcaaggcgcaatgctgcaccccccatcagctcgttcagccgttcagacacagactctgcctccatcagcagctctgcacacttctcgcgcgcctcgcagctctgcctacagcgcagattgcccttcctccacccaacagctgaacacagatttaccacgcgctgactctataa