



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

KRTAP10-3-keratin associated protein 10-3 Gene

Reference: PTXBC133677

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

KAP10.3; KAP18-3; KAP18.3; KRTAP10.3; KRTAP18-3; KRTAP18.3; keratin-associated protein 10-3; high sulfur keratin-associated protein 10.3; keratin-associated protein 18-3; keratin-associated protein 18.3; keratin associated protein 10-3

For research use only.

GENE DESCRIPTION	keratin associated protein 10-3
GENE ID	386682
ORIGIN SPECIES	Human
ACCESSIONS	BC133677
VECTOR	pDONR223
NCBI SYMBOL	KRTAP10-3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr