

HIST1H2AH-histone cluster 1, H2ah Gene

Reference: PTXBC134366

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	H2A/S; H2AFALii; H2AH; dJ86C11.1; histone H2A type 1-H; H2A histone family member; histone H2A/s; histone cluster 1, H2ah; histone cluster 1 H2A family member h		

For research use only.

GENE DESCRIPTION	histone cluster 1, H2ah
GENE ID	85235
ORIGIN SPECIES	Human
ACCESSIONS	BC134366
VECTOR	pDONR223
NCBI SYMBOL	HIST1H2AH
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctggcgctgtgaagcaggcggttaagctcgcgccaaggccaagaccgctcttctcgggctgggttcagttcccggtggccgagtgacccgctgtccgcaagggttaattatccgagcggttgagccggcgccagtgctacgtgctgctggagtagctga
ccgctgagatcctggagctggctggcaatgcggccgcgacaacaagaagaccgtatcatcccgctcactccaactggccatccgcaacgacgaggagctcaacaagctgctgggcaagtcaccatcgcgagggtgtgtcttgcctaataatcaggccgtgctgctgctaa
gaagactgagagccaccataaggccaataa