

TMEM56-transmembrane protein 56 Gene

Reference: PTXBC139737

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	transmembrane protein 56		
For research use only.			

GENE DESCRIPTION	transmembrane protein 56
GENE ID	148534
ORIGIN SPECIES	Human
ACCESSIONS	BC139737
VECTOR	pDONR223
NCBI SYMBOL	TMEM56
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagatcaacacaaactgtcatcagtggtacgtgtatcagctttttacgttttctactttgttaagttactggtttcagcaaaagtgtctcagggttcaatagctcagctcaaaaagaagattgaatggaactcaagggtagatccacatgccattcttgggtggtgtatttt
ggcctgtacattttcttattcgtatgaggtactaaagctgacccatttgggtggtcctacacttgcaaacgtgaatattgctattgctcaggctacctcatttctgatttgcattataatttgaatggaaagtgtggtgacaaattttataatgcattgtgctccctgtatgcatactacc
ttgtactgaaaaatggagtgctggcatacattgggaatttgcctgcttgcagagctttccagccggttggtaacagcgggtggttcttgaagctctgaagtatcccaagttttctaaagctatcgttatcaatggaatactcatgacagtagtattcttcatcgtgctgattgcctcaatgcttcctc
attatggcttcatgtattcgtgtatggaacagaacctacataaaggcttgaggttttaaccagttatcctgggtcattagttgtgtgttttggatgtgatgaatgtcatgtggatgatcaaaatttcaaaagggttgcataaagtcattctcacatcagacaagagaagccaaaatagcttcag
aatggaaaacttgattaa