

TRDN-triadin Gene

Reference: PTXBC139910

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CPVT5; TDN; TRISK		
For research use only.			

GENE DESCRIPTION	triadin
GENE ID	10345
ORIGIN SPECIES	Human
ACCESSIONS	BC139910
VECTOR	pDONR223
NCBI SYMBOL	TRDN
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgactgagatcactgctgaaggaaatgcactacaaccacaactgtgatagacagcaaaatggatctgtgcccaatccccgaaaagtgtgaaggagacagtcacagaagacatagtacagcgttcagctccctgcagcctggcttctgttcattgccctgataatcacgtggtc
agctgttgcctatggttattgatttagtgattacaaaactttcagcaagctctattgccagattggctcagatcctttaaactggtagctgctatggaggaaaccacggactggatctatggcttcttcttctgttatctgacatcatctcatctgaagatgaagaagatgatggtga
cgaagatagtataaaggagaatagatgagcctccttgagaaaaaagaatacacagaataagactgaaaaacaagagaacctgaaggaatacaactaaagtacacacaaagaaaaagaaaaagaaaaagtaagagaaaaagaaaaaacctgaaagaaag
caactcacaaggaaaaattgagaaaaaagaaaaccagaaacaaagacagtgaggcgaagaacagaagaagcctaagactgcagaaaagagtgaaagaaagactaaaaaggaagtgaagtggaagtggaagacaggagaaagtgaagcaaacagctgcaaaagtaaaagagtagacag
aaaacacatcaaaaccaaagaaaaggagacaaagagaagcagctgtgtcaagcatgaacagaaggacaaagccagccattccacctccctaccgacagaacaagctccagaccactccggcatcacctgccttgaaggcaagtactttttttcttaa