

Human HSP110de9 Recombinant Protein

PX-P1112-10

OVERVIEW

SIZE	10 ug
ORIGIN SPECIES	Human
FRAGMENT	Partial
PROTEIN DELIVERED WITH TAG	Yes
MOLECULAR WEIGHT WITH TAG IF ANY	43,46 kDa
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	E.coli
PURITY	>95%
PROTEIN ACCESSION	/*
FORM	liquid
BUFFER	PBS, Urea 6M
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)

PROTEIN SEQUENCE

MAHNHRHKHKLMSVVGLDVGSQSCYIAVARAGGIETIANEFSDRCTPSVISFGSKNRTIGVAAKNQQITHANNTVSNFKR
FHGRAFNDFPIQKEKENLSYDLVPLKNGGVGIKVMYMGEEHLFSVEQITAMLLTKLKETAENSLKKPVTDCVISVPSFFT
DAERRSVLDAAQIVGLNCLRLMNDMTAVALNYGIYKQDLPSLDEKPRIVVFVDMGHSAFQVSACAFNKGLKVLGTAFDP
FLGGKNFDEKLVEHFCAEFKTKYKLDAKSKIRALLRLYQECEKLKLMSSNSTDLPLNIECFMNDKDVSGKMNRSQFEEL
CAELLQKIEVPLYSLEQTHLKVEDVSAVEIVGGATRIPAVKERIAKFFGKDISTTLNADEAVARGCALQCS

More info

SPECIFICATION

GENE ID
SWISSPROTID
NCBI GENE ALIASES
UNIPROT ID
UNIPROT LINK

SYNONYMS

HSP110de9, heat shock
105kDa protein mutant de9