

Vesicular stomatitis virus VSVG IND1 Recombinant Protein

PX-P2100-100

OVERVIEW

SIZE	100 ug
ORIGIN SPECIES	Vesicular stomatitis virus
FRAGMENT	
PROTEIN DELIVERED WITH TAG	Yes
MOLECULAR WEIGHT WITH TAG IF ANY	58.44 kDa
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Insect
PURITY	25%-30%
PROTEIN ACCESSION	ACF37250.1
FORM	Frozen
BUFFER	Tris-HC 50mMl, NaCl 150mM, pH7.5 4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
STABILITY & STORAGE	

PROTEIN SEQUENCE

MKCLLYLAFLFIGVNCKFTIVFPHNQKGNWKNVPSNYHYCPSSDLNWHNDLVGTALQVKMPKSHKAIQADGWMCHASKWVTTCDFRWYGPKYITHSIRSFTP
SVEQCKESIEQTKQGTWLNPGFPPQSCGYATVTDAAIAIVQVTPHHVLVDEYTGGEWVDSQFINGKCSNDICPTVHNSTTWHSDYKVKGLCDSNLISMDITFFSED
GELSSLGKKGTGFRSNFYAYETGDKACKMQYCKHWGVRLPSGVWFEM

More info

SPECIFICATION

GENE ID

SWISSPROTID

B5AGV3

NCBI GENE ALIASES

UNIPROT ID

B5AGV3

UNIPROT LINK

<http://www.uniprot.org/uniprot/B5AGV3>

SYNONYMS

VSVG IND1, G protein ou
Glycoprotein