

Staphylococcus 6His-Protein A

PX-P2101-10mg

OVERVIEW

SIZE	10 mg
ORIGIN SPECIES	Staphylococcus
FRAGMENT	
PROTEIN DELIVERED WITH TAG	Yes
MOLECULAR WEIGHT WITH TAG IF ANY	33.83 KDa
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	E.coli
PURITY	90%
PROTEIN ACCESSION	
FORM	Frozen
BUFFER	PBS,pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)

PROTEIN SEQUENCE

MGAQHDEAAQQNAFYQVLNMPNLNADQRNGFIQSLKDDPSQSANVLGEAQLNDSQAPKADAQQNKFNKDQQSAFYEILNMPNLNEEQRNGFIQSLKDDPSQST
NVLGEAKKLNESQAPKADNNFNKEQQNAFYEILNMPNLNEEQRNGFIQSLKDDPSQSANLLAEAKKLNESQAPKADNKFNKEQQNAFYEILHLPNLNEEQRNGFI
QSLKDDPSQSANLLAEAKKLNDAQAPKADNKFNKEQQNAFYEILHLP

More info

SPECIFICATION

GENE ID

SWISSPROTID

P38507

NCBI GENE ALIASES

UNIPROT ID

P38507

UNIPROT LINK

<http://www.uniprot.org/uni/prot/P38507>

SYNONYMS

6His-Protein A