

Chinese sturgeon Follicle-stimulating hormone(Follicle) Recombinant Protein

PX-P3016-10

OVERVIEW

SIZE	10 ug
ORIGIN SPECIES	Chinese sturgeon
FRAGMENT	
PROTEIN DELIVERED WITH TAG	Yes
MOLECULAR WEIGHT WITH TAG IF ANY	39.85 kDa
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	mammalian
PURITY	90%
PROTEIN ACCESSION	
FORM	Frozen
BUFFER	PBS, pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)

PROTEIN SEQUENCE

MASVLFVLLCWAAGQCHASCALENITIGIEKDGCNCVSVNTTSCAGRCLTQADVYKSSISLYTQLVCTFKEISYVTVQLPNCPEHVDPFYTPVALSCECGQC
ATDYTDCGTLGLGSDCFSQEDDDDDKDKTHTCPPCPAPPELLGGPSVFLFPPKPKDQLMISRTPEVTCVVVDVSHEDPEVKFNWYVDGVEVHNAKTKPREEQYNS
TYRVSVLTVLHQDWLNGKEYKCKVSNKALPAPIEKTISKAKGQPREPQVYTLPPSREEMTKNQVSLTCLVKGFYPSDIAVEWESNGQPENNYKTTTPVLDSDGS
FFLYSKLTVDKSRWQQGNVFSCSVLHEALHNHYTQKSLSLSPGK

More info

SPECIFICATION

GENE ID
SWISSPROTID
NCBI GENE ALIASES
UNIPROT ID
B2CKR3
UNIPROT LINK
<http://www.uniprot.org/uni-prot/B2CKR3>

SYNONYMS

Follicle Stimulating
Hormone; FSH