

Human ChromograninA Partial [19-457] Recombinant Protein

PX-P4000-10

OVERVIEW

SIZE	10 ug
ORIGIN SPECIES	Human
FRAGMENT	Partial [19-457]
PROTEIN DELIVERED WITH TAG	Yes
MOLECULAR WEIGHT WITH TAG IF ANY	452AAs; 50.26 kDa
DELIVERY CONDITION	

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	E.coli
PURITY	0,90000000000000002
PROTEIN ACCESSION	NP_001266.1
FORM	
BUFFER	
STABILITY & STORAGE	

PROTEIN SEQUENCE

MGSSHHHHHHSSGLPVNSPMNKGDTVMKCIVEVISDTLSKPSMPVVSQECFETLRGDERILSI
LRHQNLLKELQDLALQGAKERAHQQKKHSGFEDELSEVLENQSSQAEKAEVEEPSSKDVME
KREDSKEAEKSGEATDGAEPQALPEPMQESKAEGNNQAPGEEEEEEATNTHPPASLPSQKY
PGPQAEGLDSEGLSQGLVDREKGLSAEPGWQAKREEEEEEEAEAGEEAVPEEEGPTVVNLNPH
PSLGYKEIRKGESRSEALAVDGAGKPGAEEAQDPEGKGEQEHSSQKKEEEEMAVVPQGLFRG
GKSGELEQEEERLSKEWEDSKRWSKMDQLAKELTAEKRLEGQEEEDNRDSSMKLSFRARAY
GFRGPGPQLRRGWPRSSREDSLEAGLPLQVRGYPEEKKEEEGSANRRPEDQELESLSAIEAELE
KVAHQQLALRRG

More info

UNIPROT ID
UNIPROT LINK

SYNONYMS
Chromogranin-A, CgA,
Pituitary secretory protein
I, SP-I