

Human CD58/LFA3 recombinant protein

PX-P4023-10

OVERVIEW

SIZE	10 ug
ORIGIN SPECIES	Human
FRAGMENT	
PROTEIN DELIVERED WITH TAG	Yes
MOLECULAR WEIGHT WITH TAG IF ANY	25.19kDa
DELIVERY CONDITION	any condition

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	mammalian
PURITY	0,900000000000000002
PROTEIN ACCESSION	
FORM	Lyophilized
BUFFER	50 mM Tris-HCl pH 7.5, 150 mM NaCl 4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
STABILITY & STORAGE	

PROTEIN SEQUENCE

MVAGSDAGRALGVLSVCLLHCFGFISCFSSQIYGVVYGNVTFHVPSNVPLKEVLWKKQKDK
VAELENSEFRAFSSFKNRVYLDTVSGSLTIYNLTSSDEDEYEMESPNITDTMKFFLYVLESPLSP
TLTCALTNGSIEVQCMPEYYNSHRGLIMYSWDCPMEQCKRNSTSIYFKMENDLPQKIQTLSN
PLFNTTSSILTTICPSSGHSRHRGSHHHHHH

More info

UNIPROT ID

P19256

UNIPROT LINK

<http://www.uniprot.org/uni-prot/P19256>

SYNONYMS

CD58, Ag3, Surface
glycoprotein LFA-3, LFA3,
Lymphocyte Function
Associated Antigen 3