

Human Intra Acrosomal/SP-10 recombinant protein

PX-P4029-10

OVERVIEW

SIZE	10 ug
ORIGIN SPECIES	Human
FRAGMENT	
PROTEIN DELIVERED WITH TAG	Yes
MOLECULAR WEIGHT WITH TAG IF ANY	29.7kDa
DELIVERY CONDITION	any condition

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	mammalian
PURITY	0,90000000000000002
PROTEIN ACCESSION	
FORM	Lyophilized
BUFFER	PBS, pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)

PROTEIN SEQUENCE

MNRFLLLMSLYLLGSARGTSSQPNELSGSIDHQTSTVQQLPGEFFSLENPSDAEALYETSSGLNT
LSEHGSSSEHGSSKHTVAEHTSGEHAESEHASGEPAAATEHAEGEHTVGEQPSGEQPSGEHLSGE
QPLSELESGEQPSDEQPSGEHGSGEQPSGEQASGEQPSGEHASGEQASGAPISSTSTGTILNCYT
CAYMNDQGKCLRGEGTCITQNSQQCMLKKIFEGGKLQFMVQGCENMCPSMNLFSHGTRMQII
CCRNQSFCKIGSHHHHHH

More info

UNIPROT ID

P26436

UNIPROT LINK

<http://www.uniprot.org/uni-prot/P26436>

SYNONYMS

Acrosomal protein SP 10,
Acrosomal protein SP-10,
Acrosomal vesicle protein
1, ACRV1,
ASPX_HUMAN,
D11S4365, Msa63, SP
10, Sp10, SPACA2,
Sperm protein 10