

Human Hepatitis A virus cellular receptor 2 recombinant protein

PX-P4031-10

OVERVIEW

SIZE	10 ug
ORIGIN SPECIES	Human
FRAGMENT	
PROTEIN DELIVERED WITH TAG	Yes
MOLECULAR WEIGHT WITH TAG IF ANY	23.18kDa
DELIVERY CONDITION	any condition

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	mammalian
PURITY	0,800000000000000004
PROTEIN ACCESSION	
FORM	Lyophilized
BUFFER	50 mM Tris-HCl pH 8, 150 mM NaCl 4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
STABILITY & STORAGE	

PROTEIN SEQUENCE

MFSLPFDVCVLLLLLLLLLRSSEVEYRAEVGQNAYLPCFYTPAAPGNLVPVCWGKGACPVFEC
GNVVLRTDERDVNYWTSRYWLNQDFRKGDVSLTIENVTLADSGIYCCRIQIPGIMNDEKFNLK
LVIKPAKVTPAPTLQRDFTAAAFPRMLTTRGHGPAETQTLGSLPDINLTQISTLANELRDSRLAN
DLRDSGATIRGSHHHHHH

More info

UNIPROT ID

Q8TDQ0

UNIPROT LINK

<http://www.uniprot.org/uni-prot/Q8TDQ0>

SYNONYMS

KIM-3, TIM3, TIMD3,
Tim-3, T-cell membrane
protein 3, T-cell
immunoglobulin mucin
receptor 3, T-cell
immunoglobulin and mucin
domain-containing protein
3