

Human B7RP1/ICOS/CD275 ligand recombinant protein

PX-P4032-10

OVERVIEW

SIZE	10 ug
ORIGIN SPECIES	Human
FRAGMENT	
PROTEIN DELIVERED WITH TAG	Yes
MOLECULAR WEIGHT WITH TAG IF ANY	29.56kDa
DELIVERY CONDITION	any condition

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	mammalian
PURITY	0,90000000000000002
PROTEIN ACCESSION	
FORM	Lyophilized
BUFFER	50 mM Tris-HCl pH 8, 150 mM NaCl 4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
STABILITY & STORAGE	

PROTEIN SEQUENCE

MRLGSPGLLFLLFSSLRADTQEKEVRAMVGSDVELSCACPEGSRFDLNDVYVYWQTSESKTV
VTYHIPQNSSLENVDSRYRNRALMSPAGMLRGDFSLRLFNVTTPQDEQKFHCLVLSQSLGFQEV
LSVEVTLHVAANFSVPVVSAPHSPSQDELFTCTTSINGYPRPNVYWINKTDNSLLDQALQNDT
VFLNMRGLYDVVSVLRIARTPSVNIGCCCIENVLLQQNLTVGSQTGNDIGERDKITENPVSTGEK
NAATWSGSHHHHHH

More info

UNIPROT ID

O75144

UNIPROT LINK

<http://www.uniprot.org/uni-prot/O75144>

SYNONYMS

CD275, B7-H2, B7H2,
B7RP-1, B7RP1, GL50,
ICOS-L, ICOSL, LICOS,
B7-like protein Gl50, B7-
related protein 1