

Human CD54/ICAM1 recombinant protein

PX-P4033-10

OVERVIEW

SIZE	10 ug
ORIGIN SPECIES	Human
FRAGMENT	
PROTEIN DELIVERED WITH TAG	Yes
MOLECULAR WEIGHT WITH TAG IF ANY	25.62kDa
DELIVERY CONDITION	any condition

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	mammalian
PURITY	0,8000000000000004
PROTEIN ACCESSION	
FORM	Lyophilized
BUFFER	50 mM Tris-HCl pH 8, 150 mM NaCl 4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
STABILITY & STORAGE	

PROTEIN SEQUENCE

MAPSSRPALPALLVLLGALFPGPGNAQTSVSPSKVILPRGGSVLVTCSTSCDQPKLLGIETPLPKKELLPG
NNRKVYELSNVQEDSQPMCYSNCPDGQSTAKTFLTVYWTPERVELAPLPSWQPVGKNLTLRCQVEGGAP
RANLTVVLLRGEKELKREPAVGEPAEVTTTTLVRRDHHGANFSCRTDLRPQGLELFENTSAPYQLQTFV
LPATPPQLVSPRVLEVDQTQGTVVCSLDGLFPVSEAQVHLALGDQRLNPTVTYGNDSFSAKASVSVTAEDG
TQRLTCAVILGNQSQETLQTVTIYSFPAPNVILTKPEVSEGTEVTVKCEAHPRAKVTLNGVPAQPLGPRAQL
LLKATPEDNGRSFSCSATLEVAGQLIHKNTRELRLVLYGPRLDERDCPGNWTWPENSQQTPMCQAWGNPL
PELKCLKDGTFPLPIGESVTVTRDLEGTYLCRARSTQGEVTRKVTNVLSPRYEGSHHHHHH

More info

UNIPROT ID

P05362

UNIPROT LINK

<http://www.uniprot.org/uni-prot/P05362>

SYNONYMS

CD54, BB2, P3.58, Human
Rhinovirus Receptor,
Major group rhinovirus
receptor, Intercellular
Adhesion Molecule 1

For research use only.