

Human CD102/ICAM2 recombinant protein

PX-P4034-10

OVERVIEW

SIZE	10 ug
ORIGIN SPECIES	Human
FRAGMENT	
PROTEIN DELIVERED WITH TAG	Yes
MOLECULAR WEIGHT WITH TAG IF ANY	25.62kDa
DELIVERY CONDITION	any condition

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	mammalian
PURITY	0,800000000000000004
PROTEIN ACCESSION	
FORM	Lyophilized
BUFFER	50 mM Tris-HCl pH 8, 150 mM NaCl 4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
STABILITY & STORAGE	

PROTEIN SEQUENCE

MSSFGYRTLTVLFTLICCPGSDEKVFEVHVRPKKLAVEPKGSLEVNCSTTCNQPEVGGLETSL
DKILLDEQAQWKHYLVSNISHDTVLQCHFTCSGKQESMNSNVSVYQPPRQVILTLQPTLVAVG
KSFTIECRVPTVEPLDSLTLFLFRGNETLHYETFGKAAPAPQEATATFNSTADREDGHRNFSCL
AVLDLMSRGGNIFHKHSAPKMLEIYEPVSDSQGSHHHHHH

More info

UNIPROT ID

P13598

UNIPROT LINK

<http://www.uniprot.org/uni-prot/P13598>

SYNONYMS

ICAM2, CD102, CD102
antigen, ICAM 2,
Intercellular adhesion
molecule 2, Ly60