

Human T-cell immunoreceptor with Ig and ITIM domains (TIGIT) recombinant protein

PX-P4036-10

OVERVIEW

SIZE	10 ug
ORIGIN SPECIES	Human
FRAGMENT	
PROTEIN DELIVERED WITH TAG	Yes
MOLECULAR WEIGHT WITH TAG IF ANY	16.36kDa
DELIVERY CONDITION	any condition

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	mammalian
PURITY	0,400000000000000002
PROTEIN ACCESSION	
FORM	Lyophilized
BUFFER	50 mM Tris-HCl pH 8, 150 mM NaCl 4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
STABILITY & STORAGE	

PROTEIN SEQUENCE

MRWCLLLIWAQGLRQAPLASGMMTGTIETTGNI~~SA~~EKGGSII~~L~~QCHLSSTTAQVTQVNWEQQ
DQLLAICNADLGWHISPSFKDRVAPGPG~~L~~GLTLQSLTVNDTGEYFCIYHTYPDGTYTGRIFLEV
LESSVAEHGARFQIPGSHHHHHH

More info

UNIPROT ID

Q495A1

UNIPROT LINK

<http://www.uniprot.org/uniprot/Q495A1>

SYNONYMS

TIGIT, VSIG9, VSTM3,
V-Set And
Immunoglobulin Domain
Containing 9, V-Set And
Transmembrane Domain
Containing 3