

Human Tumor necrosis factor receptor superfamily member 1A (TNFR1A) recombinant protein

PX-P4037-10

OVERVIEW

SIZE	10 ug
ORIGIN SPECIES	Human
FRAGMENT	
PROTEIN DELIVERED WITH TAG	Yes
MOLECULAR WEIGHT WITH TAG IF ANY	24.35kDa
DELIVERY CONDITION	any condition

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	mammalian
PURITY	0,90000000000000002
PROTEIN ACCESSION	
FORM	Lyophilized
BUFFER	50 mM Tris-HCl pH 8, 150 mM NaCl 4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
STABILITY & STORAGE	

PROTEIN SEQUENCE

MGLSTVPDLLLPLVLELLVGIYPSGVIGLVPHLGDREKRDSVCPQGKYIHPQNNSICCTKCHK
GTYLYNDPCPGPGQDTCRECESGSFTASENHLRHCLSCSKCRKEMGQVEISSCTVDRDTVCGC
RKNQYRHYWSENLFQCFNCSLCLNGTVHLSCQEKQNTVCTCHAGFFLRENECVSCSNCKKSL
ECTKLCLPQIENVKGTEDSGTTGSHHHHHH

More info

UNIPROT ID

P19438

UNIPROT LINK

<http://www.uniprot.org/uniprot/P19438>

SYNONYMS

CD120A, P55, TBP1, FPF,
TNF-R, TNF-R-I, TNF-
R55, TNFAR, TNFR1,
TNFR55, TNFR60, P55-R,
P60, Tumor necrosis factor
receptor 1, Tumor necrosis
factor-binding protein 1