

Human CD270/TNFRSF14/HVEM recombinant protein

PX-P4039-10

OVERVIEW

SIZE	10 ug
ORIGIN SPECIES	Human
FRAGMENT	
PROTEIN DELIVERED WITH TAG	Yes
MOLECULAR WEIGHT WITH TAG IF ANY	22.46kDa
DELIVERY CONDITION	any condition

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	mammalian
PURITY	0,40000000000000002
PROTEIN ACCESSION	
FORM	Lyophilized
BUFFER	PBS, pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)

PROTEIN SEQUENCE

MEPPGDWGWPPWRSTPKTDVLRLLVLYLTFLGAPCYAPALPSCKEDEYPVGSECCPKCSPGYRV
KEACGELTGTVCEPCPPGTYYIAHLNGLSKCLQCQMCDPAMGLRASRNCSTENAVCGCSPGH
FCIVQDGDHCAACRAYATSSPGQVRVQKGGTESQDTLCQNCPPGTFSPNGTLEECQHQTCSW
LVTKAGAGTSSSHWVGSHHHHHH

More info

UNIPROT ID

Q92956

UNIPROT LINK

<http://www.uniprot.org/uni-prot/Q92956>

SYNONYMS

CD270, TR2, ATAR,
HVEA, HVEM, LIGHTR,
Herpesvirus Entry
Mediator, Herpes virus
entry mediator A, Tumor
necrosis factor receptor-like
2