

Human CD80/B7-1 recombinant protein

PX-P4040-10

OVERVIEW

SIZE	10 ug
ORIGIN SPECIES	Human
FRAGMENT	
PROTEIN DELIVERED WITH TAG	Yes
MOLECULAR WEIGHT WITH TAG IF ANY	28.52kDa
DELIVERY CONDITION	any condition

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	mammalian
PURITY	0,9000000000000000
PROTEIN ACCESSION	
FORM	Lyophilized
BUFFER	PBS, pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)

PROTEIN SEQUENCE

MGHTRRQGTSPSKCPYLNFFQLVLGLSHFCSGVIHVTKEVKEVATLSCGHNVSVLELAQTR
IYWQKEKKMVLTMMSGDMNIWPEYKNRTIFDITNNLSIVILALRPSDEGTYECVVLKYEKDAF
KREHLAEVTLVKADFPTPSISDFEIPTSNIRRIICSTSGGFPEPHLSWLENGEELNAINTTVSQDP
ETELYAVSSKLDNFMTTNHSFMCLIKYGHLRVNQTFNWNTTKQEHFPDNGSHHHHHH

More info

UNIPROT ID

P33681

UNIPROT LINK

<http://www.uniprot.org/uni-prot/P33681>

SYNONYMS

CD80, CD28LG,
CD28LG1, BB1, B7-1
Antigen, T-lymphocyte
activation antigen CD80,
CTLA-4 counter-receptor
B7.1