

Human CD28 recombinant protein

PX-P4044-10

OVERVIEW

SIZE	10 ug
ORIGIN SPECIES	Human
FRAGMENT	
PROTEIN DELIVERED WITH TAG	Yes
MOLECULAR WEIGHT WITH TAG IF ANY	42.63kDa
DELIVERY CONDITION	any condition

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	mammalian
PURITY	0,69999999999999996
PROTEIN ACCESSION	
FORM	Lyophilized
BUFFER	PBS, pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)

PROTEIN SEQUENCE

MLRLLLALNLFPSIQVTGNKILVKQSPMLVAYDNAVNLSCKYSYNLFSREFRASLHKGLDSAVEV
CVVYGNYSQQQLQVYSKTGFNCDGKLGNESVTFYLQNLVYNQTDIYFCKIEVMYPPPYLDNEKSN
GTIIHVKGKHLCPSPFLPGPSKPDKTHTCPAPPELLGGPSVFLFPPKPKDQLMISRTPEVTCVVVD
VSHEDPEVKFNWYVDGVEVHNAKTKPREEQYNSTYRVVSVLTVLHQDWLNGKEYKCKVSNKA
LPAPIEKTISKAKGQPREPQVYTLPPSREEMTKNQVSLTCLVKGFYPSDIAVEWESNGQPENNYKT
TPPVLDSDGSFFLYSKLTVDKSRWQQGNVFSFCSVLHEALHNHYTQKSLSLSPGK

More info

UNIPROT ID

P10747

UNIPROT LINK

<http://www.uniprot.org/uniprot/P10747>

SYNONYMS

CD28, T-cell-specific
surface glycoprotein CD28,
Tp44, T-Cell-Specific
Surface Glycoprotein,
MGC138290