

Human CD126/IL6R/Interleukin-6 receptor subunit alpha recombinant protein

PX-P4046-10

OVERVIEW

SIZE	10 ug
ORIGIN SPECIES	Human
FRAGMENT	
PROTEIN DELIVERED WITH TAG	Yes
MOLECULAR WEIGHT WITH TAG IF ANY	41.12kDa
DELIVERY CONDITION	any condition

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	mammalian
PURITY	0,800000000000000004
PROTEIN ACCESSION	
FORM	Lyophilized
BUFFER	PBS, pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)

PROTEIN SEQUENCE

MLAVGCALLAALLAAPGAALAPRRCPAQEVARGVLTSLPGDSVTLTCPGVEPEDNATVHWVLRKPAAGS
HPSRWAGMGRLLRLRSVQLHDSGNYSCYRAGRPAQTVHLLVDVPPEEPQLSCFRKSPLSNVVCEWGRST
PSLTTKAVLLVRKFQNSPAEDFQEPQYSQESQKFSCQLAVPEGDSSFYIVSMCVASSVGSKFSTQTFQGC
GILQPDPPANITVTAVARNPRWLSVTWQDPHSWNSSFYRLRFELRYRAERSKFTTTWMVKDLQHHCVIHD
AWSGLRHVVQLRAQEEFGQGEWSEWSPEAMGTPWTESRPPAENEVSTPMQALTTNKDDDDNILFRDSAN
ATSLPVQDSSSVPLPGSHHHHHH

More info

UNIPROT ID

P08887

UNIPROT LINK

<http://www.uniprot.org/uniprot/P08887>

SYNONYMS

CD126, gp80, IL6R, IL-6R-1, IL-6R-Alpha, IL6RA, IL-6 receptor subunit alpha, Membrane glycoprotein 80

For research use only.