

Human BNP peptide/Natriuretic peptides B

PX-P4047-10

OVERVIEW

SIZE	10 ug
ORIGIN SPECIES	Human
FRAGMENT	
PROTEIN DELIVERED WITH TAG	Yes
MOLECULAR WEIGHT WITH TAG IF ANY	28.85kDa
DELIVERY CONDITION	any condition

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	E.coli
PURITY	0,8499999999999998
PROTEIN ACCESSION	
FORM	Lyophilized
BUFFER	50 mM Tris-HCl pH 8, 150 mM NaCl 4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
STABILITY & STORAGE	

PROTEIN SEQUENCE

MGSHHHHHHSGHPLGSPGSASDLETSGLQEQRNHLQGKLSSELQVEQTSLEPLQESPRPTGVWK
SREVATEGIRGHRKMLVLYTLRAPR

More info

UNIPROT ID

P16860

UNIPROT LINK

<http://www.uniprot.org/uni-prot/P16860>

SYNONYMS

NT-PROBNP, BNP,
Natriuretic peptide B, GC-
B, B-Type Natriuretic
Peptide, Ventricular
Natriuretic Peptide,
Gamma-brain natriuretic
peptide, Brain natriuretic
peptide 32