

Human Golgi membrane protein 1 (GP73) recombinant protein

PX-P4051-10

OVERVIEW

SIZE	10 ug
ORIGIN SPECIES	Human
FRAGMENT	
PROTEIN DELIVERED WITH TAG	Yes
MOLECULAR WEIGHT WITH TAG IF ANY	42.35kDa
DELIVERY CONDITION	any condition

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	E.coli
PURITY	0,800000000000000004
PROTEIN ACCESSION	
FORM	Lyophilized
BUFFER	50 mM Tris-HCl pH 8, 150 mM NaCl 4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
STABILITY & STORAGE	

PROTEIN SEQUENCE

MGSHHHHHHSGVDLQTRIMELEGRVRRAAAERGAVELKKNEFQGELEKQREQLDKIQSSHNFQL
ESVNKLYQDEKAVLVNNITTGERLIRVLQDQLKTLQRNYGRLQQDVLQFQKNQTNLERKFSYDLS
QCINQMKEVKEQCEERIEEVTKKGNEAVASRDLENNDQRQLQALSEPQRLQAAGLPHTTEVPQ
GKGNVLGNSKSQTPAPSSEVVLDKSRQVEKEETNEIQVVNEEPQRDRLPQEPGREQVVEDRPVGG
RGFGGAGELGQTPQVQAALSVSQENPEMEGPERDQLVIPDGQEEEQEAAGEGRNQKLRGEDDY
NMDENAESETDKQAALAGNDRNIDVFNVEDQKRDITINLLDQREKRNHTL

More info

UNIPROT ID

Q8NBJ4

UNIPROT LINK

<http://www.uniprot.org/uni/prot/Q8NBJ4>

SYNONYMS

GOLM1, GOLPH2,
C9orf155, Golgi Membrane
Protein 1, Golgi
Phosphoprotein 2, GP73