

Rhesus monkey IL37/Interleukin 1 Zeta (IL1z) recombinant protein

PX-P4063-10

OVERVIEW

SIZE	10 ug
ORIGIN SPECIES	Rhesus monkey
FRAGMENT	Full-length
PROTEIN DELIVERED WITH TAG	Yes
MOLECULAR WEIGHT WITH TAG IF ANY	27.71kDa
DELIVERY CONDITION	

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	E.coli
PURITY	0,9499999999999996
PROTEIN ACCESSION	
FORM	
BUFFER	
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)

PROTEIN SEQUENCE

MGSHHHHHHSGLVPRGSMNNSTLKMSFVGNSGVKTGSEDWEKDEPQCYSEKDEPQCYSE
DPAGSPLEPGPSLPSMNFVHTSPKVKNLNPKKFSIHDQDHKVLVVDSGNLIAVPDKNYIRPEIFF
ALASSLSASAEGSPILLGVSKGEFCLSCDKDKGQSHPSLQLKKKKLMKLAALKESARRPFIF
YRAQVGSRTLESAAHPGWFICTSCNENEPVGVTDKFENRKHIEFSFQPVCKAEMSPSEVSD

More info

UNIPROT ID

UNIPROT LINK

SYNONYMS

IL1F7, IL37, FIL1,
FIL1(ZETA), FIL1Z, IL-
1F7, IL-1H4, IL-1RP1,
IL1H4, IL1RP1,
Interleukin 37, Interleukin-
1-Related Protein,
Interleukin 1 Family,
Member 7

For research use only.