

Mouse Monoclonal antibody to human 9G8 (Splicing factor)

OVERVIEW

CATALOG NO.	PTGX-9G8-1C6
SIZE	100µL
CLONALITY	Monoclonal Antibody
HOST	Mouse
CLONE NAME	2G8-1C6
ISOTYPE	IgG1, k

APPLICATIONS AND STORAGE

REACTIVITY	Human, Mouse
APPLICATIONS	E, WB, IC IP
FORM	liquid/frozen
DELIVERY CONDITION	Dry Ice
STABILITY & STORAGE	Store at -20°C. Aliquot to avoid repeated freezing and thawing.