

Mouse Monoclonal antibody to human Actin

OVERVIEW

CATALOG NO.	PTGX-ACT-2D7
SIZE	100µL
CLONALITY	Monoclonal Antibody
HOST	Mouse
CLONE NAME	1ACT-2D7
ISOTYPE	IgG2a, k

APPLICATIONS AND STORAGE

REACTIVITY	Human, Mouse, Rat
APPLICATIONS	E, WB, IC IP
FORM	liquid/frozen
DELIVERY CONDITION	Dry Ice
STABILITY & STORAGE	Store at -20°C. Aliquot to avoid repeated freezing and thawing.