

Mouse Monoclonal antibody to human Cyclin L1

OVERVIEW

CATALOG NO.	PTGX-CYL-2F10
SIZE	100µL
CLONALITY	Monoclonal Antibody
HOST	Mouse
CLONE NAME	3CY-2F10
ISOTYPE	IgG2a, k

APPLICATIONS AND STORAGE

REACTIVITY	Human, Mouse, Rat
APPLICATIONS	E, WB, IC, IHC
FORM	liquid/frozen
DELIVERY CONDITION	Dry Ice
STABILITY & STORAGE	Store at -20°C. Aliquot to avoid repeated freezing and thawing.