

Mouse Monoclonal antibody to human DAX-1(NROB1)

OVERVIEW

CATALOG NO.	PTGX-DAX-2F4
SIZE	100µL
CLONALITY	Monoclonal Antibody
HOST	Mouse
CLONE NAME	1DA-2F4
ISOTYPE	IgG1, k

APPLICATIONS AND STORAGE

REACTIVITY	Human, Mouse, Rat
APPLICATIONS	E, WB, IC, IHC, IP
FORM	liquid/frozen
DELIVERY CONDITION	Dry Ice
STABILITY & STORAGE	Store at -20°C. Aliquot to avoid repeated freezing and thawing.