



ProteoGenix SAS
Espace Européen de l'Entreprise
15 rue de La Haye
67300 Schiltigheim, France
Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88
contact@proteogenix.fr
www.ProteoGenix.science

Mouse Monoclonal antibody to human Histone H3 Dimethyl (K)-Phospho (S),

OVERVIEW

| | |
|--------------------|---------------------|
| CATALOG NO. | PTGX-H3-2H12 |
| SIZE | 100µL |
| CLONALITY | Monoclonal Antibody |
| HOST | Mouse |
| CLONE NAME | 51TA-2H12 |
| ISOTYPE | IgG1, k |

APPLICATIONS AND STORAGE

| | |
|--------------------------------|---|
| REACTIVITY | Human, Mouse, Rat |
| APPLICATIONS | E, WB, IC, IP |
| FORM | liquid/frozen |
| DELIVERY CONDITION | Dry Ice |
| STABILITY & STORAGE | Store at -20°C. Aliquot to avoid repeated freezing and thawing. |

For research use only.