

Mouse Monoclonal antibody to human Histone H3 (Phospho-Ser10)

OVERVIEW

CATALOG NO.	PTGX-H3-4C8
SIZE	100µL
CLONALITY	Monoclonal Antibody
HOST	Mouse
CLONE NAME	3HH3-4C8
ISOTYPE	IgG2a, k

APPLICATIONS AND STORAGE

REACTIVITY	Human, Mouse, Rat
APPLICATIONS	E, WB, IC, IP
FORM	liquid/frozen
DELIVERY CONDITION	Dry Ice
STABILITY & STORAGE	Store at -20°C. Aliquot to avoid repeated freezing and thawing.