

Mouse Monoclonal antibody to human TBP

OVERVIEW

CATALOG NO.	PTGX-TBP-2C1
SIZE	100µL
CLONALITY	Monoclonal Antibody
HOST	Mouse
CLONE NAME	5TF1-2C1
ISOTYPE	IgG2b, k

APPLICATIONS AND STORAGE

REACTIVITY	Human, Mouse, Chicken
APPLICATIONS	E, WB, IC, IHC, IP
FORM	liquid/frozen
DELIVERY CONDITION	Dry Ice
STABILITY & STORAGE	Store at -20°C. Aliquot to avoid repeated freezing and thawing.