

## Mouse Monoclonal antibody to human TBP

### OVERVIEW

---

<b>CATALOG NO.</b>	PTGX-TBP-3G3
<b>SIZE</b>	100µL
<b>CLONALITY</b>	Monoclonal Antibody
<b>HOST</b>	Mouse
<b>CLONE NAME</b>	3TF1-3G3
<b>ISOTYPE</b>	IgG1, k

### APPLICATIONS AND STORAGE

---

<b>REACTIVITY</b>	Human, Mouse, Xenopus, Chicken, Zebrafish
<b>APPLICATIONS</b>	E, WB, IC, IP
<b>FORM</b>	liquid/frozen
<b>DELIVERY CONDITION</b>	Dry Ice
<b>STABILITY &amp; STORAGE</b>	Store at -20°C. Aliquot to avoid repeated freezing and thawing.