

Cleaved-Caspase-6 p18 (D162) Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX10002	SYNONYMS
SIZE	100 ug	CASP6; MCH2; Caspase-6; CASP-6; Apoptotic protease
HOST	Rabbit	Mch-2

PRODUCT INFORMATION

GENE NAME	CASP6
PROTEIN NAME	Caspase-6
IMMUNOGEN	Synthesized peptide derived from the Internal region of mouse Caspase-6 p18. at AA range: 120-200

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Mouse,Rat
TESTED APPLICATIONS	WB,IHC,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.