

Cleaved-Caspase-6 p18 (D179) Polyclonal Antibody

OVERVIEW

| | | |
|--------------------|----------|--|
| CATALOG NO. | PTX10003 | SYNONYMS |
| SIZE | 100 ug | CASP6; MCH2; Caspase-6; CASP-6; Apoptotic protease |
| HOST | Rabbit | Mch-2 |

PRODUCT INFORMATION

| | |
|---------------------|--|
| GENE NAME | CASP6 |
| PROTEIN NAME | Caspase-6 |
| IMMUNOGEN | Synthesized peptide derived from the Internal region of human Caspase-6 p18. at AA range: 90-170 |

APPLICATIONS AND STORAGE

| | |
|--------------------------------|--|
| SPECIES REACTIVITY | Human,Rat |
| TESTED APPLICATIONS | WB,IHC,ELISA |
| RECOMMENDED DILUTION | Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications. |
| DELIVERY CONDITIONS | Blue Ice |
| STABILITY & STORAGE | PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles |

For research use only.